




GH3+

Disruptive Technology in Safe Patient Handling & Mobility



GH3+ Ceiling Lift System

Always ready for use

Continuous charge rail technology eliminates the need for a charging dock.

Follow-up Intervention

CLM Online. View the day and time of lifts, how long transfers are taking and the height and weight of each transfer anywhere at anytime via wifi connection and online dashboard.

Do more

The Class III Scale is the first approved medical scale on the market that is integrated into a ceiling lift, fulfilling the conditions for medical weighing.

Early Mobility

The Trainer Module, with active assisted therapeutic activity capabilities, makes it possible to start rehabilitation with dynamic weight relief in the patient's room.

Rehabilitation Exercises

The Positioning Lock system secures the GH3 and/or the traverse rail in a given position in the full-room covering rail system. This allows the GH3+ to be used for the rehabilitation and retraining of both bed-bound and more mobile users.

Scan me



View a video
of the GH3+
Trainer Module

guldmann.com | info@guldmann.net

Time to care