



Guldmann™



# Disruptive Technology in Safe Patient Handling & Mobility

## **GH3** + **Trainer Module**

| Time to care |

### Improve mobility and physical function with dynamic body weight support

#### Do more than lifting up and down

The Trainer Module, with active assisted therapeutic activity capabilities, makes it possible to start rehabilitation with dynamic weight relief in the patient's room. The Trainer Module adds another option to the many application possibilities that the Guldmann GH3+ ceiling hoist offers besides lifting up and down.

#### Dynamic weight relief

The GH3+ Trainer Module will dynamically unload the patient's body weight up to 220lbs. The caregiver selects the amount of weight needed to be unloaded for the patient and the Trainer module will do the rest.

[guldmann.com](http://guldmann.com) | [info@guldmann.net](mailto:info@guldmann.net)



GH3+ Trainer Module demonstration.