



Guldmann™

# Disruptive Technology in Safe Patient Handling & Mobility



## GH3+ CLM Online <sup>wifi</sup>

| Time to care |

### A tool for professionals to increase, enhance & upgrade control over their SPHM program


#### Using the data

CLM Online was developed to create an overview and distribute data of use within large facilities with multiple GH3+ ceiling lift systems in multiple buildings, wards and rooms.

#### Follow-up Intervention

In real-time and from anywhere you have wi-fi access, view the day and time of lifts, how long transfers are taking and the height and weight of each transfer. Observing lift use and trends provides opportunity for SPHM Program leaders to objectively quantify and establish lift use expectations.

Contact us today for a live demonstration .



[guldmann.com](http://guldmann.com) | [info@guldmann.net](mailto:info@guldmann.net)